Ispitni zadaci za kolegij Difrakcija u polikristalnom materijalu 22. rujna 2021.

Zadatak 1.

Izvedite jednadžbu (1) prema tekstu i slici ispod.

1. Introduction

Flat-plate asymmetric reflection geometry is a useful arrangement for diffraction studies of surface phenomena and *in situ* experimentation, and is the only geometry possible when using a flat-plate sample with a wide-angle position-sensitive detector, such as the INEL CPS120 in laboratory instruments or the Mythen detector (Schmitt *et al.*, 2003) as used at a number of synchrotron facilities. As the use of this geometry becomes more widespread, it becomes increasingly important to correctly model the effects of the instrument on the measured diffraction patterns. Toraya & Yoshino (1994) have demonstrated the corrections for a pseudo-parallel beam using a modified Thompson-Cox-Hastings model (Thompson *et al.*, 1987); however, the geometry of the instrument described by Toraya and Yoshino included a Ge(111) analyser crystal in the diffracted beam, which has a significant impact on the observed peak widths and shapes. In the configuration described here, there are no diffracted beam optics; therefore, it is the size and intensity distribution of the footprint of the incident beam on the sample that governs the observed peak widths and shapes, and hence the overall peak resolution. This paper describes the development of a fundamental parameters model which uses the known instrument settings to calculate the instrument component of width.

In this asymmetric geometry, the observed intensity has a different variation as a function of diffraction angle compared with conventional Bragg-Brentano arrangements (Toraya *et al.*, 1993), and the effect of sample displacement is greatly enhanced (Masson *et al.*, 1996). In addition to these effects, the width of the peak profiles changes significantly with increasing diffraction angle (see Fig. 1) owing to the 2θ dependence of the apparent size of the beam on the sample, which is due to the absence of any diffracted beam optics.

Assuming a parallel incident beam, the width, w , of the diffracted beam, in degrees, is given by

$$w = \frac{b}{\sin(\alpha)} \frac{\sin(2\theta - \alpha)}{R} \frac{180}{\pi}, \quad (1)$$

where b is the beam height, α is the angle of the sample with respect to the incident beam (beam angle), R is the radius of the detector and 2θ is the diffraction angle. The first term is the actual length of the beam on the sample (beam footprint) and the second term modifies the beam footprint to give the apparent length as seen at a given diffraction angle. The width of the diffracted beam goes through a maximum at an angle perpendicular to the sample surface. All


Figure 1

A schematic diagram showing the effect of diffraction angle (2θ) on the width (w) of the diffracted beam in flat-plate asymmetric reflection geometry. The height of the beam is given by b ; the angle of the sample with respect to the incident beam (beam angle) is given by α . The highlighted section of the sample indicates the beam footprint.

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